

## Quick Start Guide



by  neurolens®

Follow us on:



[facebook.com/Neurolenses](https://facebook.com/Neurolenses)



[@neurolenses](https://@neurolenses)



[@neurolens](https://@neurolens)



[youtube.com/@neurolens](https://youtube.com/@neurolens)

[neuroLens.com](https://neuroLens.com)

# N<sup>3</sup> Quick Start Guide

1



Obtain your PIN from the Practice Profile --> User Management section of the Neurolens Portal

2



Press the Power button on top of the device.

3




Locate the "Confirm" button on the side of the device.

4



Place the device on your head and adjust the dial to fit.

5



Follow the instructions to ensure device is ready for use.

6



Start changing lives!

To access the full user manual, or to request a printed copy, go to [neurolens.com/n3-manual](https://neurolens.com/n3-manual)